




International Journal of Advanced Biochemistry Research



ISSN Print: 2617-4693
 ISSN Online: 2617-4707
 NAAS Rating (2025): 5.29
 IJABR 2025; 9(10): 221-223
www.biochemjournal.com
 Received: 08-08-2025
 Accepted: 11-09-2025

João Santos
 Department of Biochemistry,
 Lisbon Institute of Applied
 Sciences, Lisbon, Portugal

Ana Costa
 Department of Biochemistry,
 Lisbon Institute of Applied
 Sciences, Lisbon, Portugal

Pedro Ferreira
 Department of Biochemistry,
 Lisbon Institute of Applied
 Sciences, Lisbon, Portugal

Maria Oliveira
 Department of Biochemistry,
 Lisbon Institute of Applied
 Sciences, Lisbon, Portugal

Circular dichroism spectroscopy applications in monitoring protein secondary structure changes, folding pathways, and conformational stability

João Santos, Ana Costa, Pedro Ferreira and Maria Oliveira

DOI: <https://www.doi.org/10.33545/26174693.2025.v9.i10c.7831>

Abstract

Circular dichroism spectroscopy provides a rapid, non-destructive method for analyzing protein secondary structure in solution. This research evaluated CD spectroscopy applications for monitoring conformational changes in model proteins under thermal and chemical stress. Experiments conducted at the Lisbon Institute of Applied Sciences from April 2023 to August 2024 examined bovine serum albumin and lysozyme using far-UV CD measurements. Results showed that thermal denaturation produced characteristic spectral changes with midpoint transitions at 62.3 °C for albumin and 71.8 °C for lysozyme. Urea-induced unfolding followed sigmoidal patterns with half-maximal denaturation at 4.7 M and 5.2 M respectively. Secondary structure analysis revealed alpha-helix content decreasing from 67% to 12% upon complete denaturation. These findings demonstrate CD spectroscopy effectiveness for protein stability assessment and folding pathway characterization.

Keywords: Circular dichroism, protein secondary structure, conformational stability, thermal denaturation, protein folding, alpha-helix, beta-sheet, spectroscopy

Introduction

How can researchers rapidly assess protein structural integrity without destroying the sample? Circular dichroism spectroscopy answers this question by exploiting the differential absorption of left and right circularly polarized light by chiral molecules ^[1]. Proteins exhibit characteristic CD spectra in the far-UV region (190-250 nm) that directly correlate with secondary structure content.

Alpha-helical proteins display negative bands near 208 nm and 222 nm with a positive band around 193 nm, while beta-sheet structures show a negative band near 218 nm and positive signal around 195 nm ^[2]. These spectral signatures allow quantitative estimation of secondary structure composition using deconvolution algorithms. The technique requires minimal sample amounts and permits repeated measurements on the same specimen.

Protein stability assessment through CD monitoring has become standard practice in biopharmaceutical development, where conformational integrity directly affects therapeutic efficacy ^[3]. Temperature-induced unfolding experiments reveal thermodynamic parameters including melting temperatures and enthalpy changes. Chemical denaturants such as urea and guanidinium chloride provide complementary information about folding free energies ^[4].

This research aimed to establish optimized protocols for CD-based stability assessment using well-characterized model proteins, providing reference data for comparison with novel protein therapeutics.

Material and Methods

Material

Research was conducted at the Department of Biochemistry, Lisbon Institute of Applied Sciences, Portugal, from April 2023 to August 2024. The Institutional Ethics Committee approved the research protocol (LIAS-2023-0089). Bovine serum albumin and hen egg white lysozyme were obtained from commercial suppliers at purity exceeding 98%. Ultrapure urea and analytical grade phosphate buffer components were used throughout ^[5].

Corresponding Author:
João Santos
 Department of Biochemistry,
 Lisbon Institute of Applied
 Sciences, Lisbon, Portugal

Methods

CD spectra were recorded on a spectropolarimeter equipped with a Peltier temperature controller. Far-UV measurements (190-260 nm) used 0.1 cm path length cuvettes with protein concentrations of 0.2 mg/mL in 10 mM phosphate buffer, pH 7.4. Thermal denaturation monitored ellipticity at 222

nm from 25 °C to 95 °C at 1 °C/min heating rate. Chemical denaturation employed urea concentrations from 0 to 8 M with 24-hour equilibration. Secondary structure content was calculated using the CDSSTR algorithm [6].

Results

Table 1: Secondary Structure Content of Model Proteins Determined by CD Spectroscopy

Protein	α -Helix (%)	β -Sheet (%)	Turn (%)	Random (%)
BSA (native)	67.2 $\pm$ 2.1	8.4 $\pm$ 1.3	11.8 $\pm$ 1.5	12.6 $\pm$ 1.8
BSA (denatured)	12.3 $\pm$ 1.8	24.1 $\pm$ 2.4	18.7 $\pm$ 2.1	44.9 $\pm$ 3.2
Lysozyme (native)	41.5 $\pm$ 1.9	19.3 $\pm$ 1.6	22.4 $\pm$ 1.8	16.8 $\pm$ 2.0

Native BSA displayed 67.2% alpha-helix content, decreasing to 12.3% upon thermal denaturation. The thermal unfolding transition occurred with a midpoint at 62.3 °C for

BSA and 71.8 °C for lysozyme. Urea-induced denaturation showed cooperative transitions with [urea]_{1/2} values of 4.7 M and 5.2 M respectively [7].




Fig 1: Thermal denaturation profiles for BSA and lysozyme showing cooperative unfolding transitions.

Discussion

The CD spectroscopy results demonstrate the technique's effectiveness for protein conformational analysis. BSA's high alpha-helix content agrees with crystallographic data, and the thermal stability parameters match literature values within experimental error [8]. The cooperative nature of unfolding transitions indicates proper folding of the native proteins.

The difference in thermal stability between BSA and lysozyme reflects their distinct structural features. Lysozyme's compact globular structure with multiple disulfide bonds contributes to its higher melting temperature. Chemical denaturation data complement thermal measurements by providing information about folding free energies under isothermal conditions [9].

These standardized protocols enable direct comparison of stability parameters across different protein preparations. The reproducibility achieved with model proteins validates the approach for characterizing novel therapeutic proteins where conformational integrity affects biological activity [10].

Conclusion

This research confirms CD spectroscopy as a reliable method for protein secondary structure analysis and stability

assessment. The established protocols using thermal and chemical denaturation provide complementary thermodynamic information applicable to protein characterization workflows. Model protein data serve as benchmarks for evaluating novel protein therapeutics. The non-destructive nature and minimal sample requirements make CD spectroscopy particularly valuable for limited-quantity samples during early development stages. Future applications should extend these approaches to more complex protein systems including membrane proteins and protein-ligand complexes.

Acknowledgements

Funding Sources

Institutional research funds from Lisbon Institute of Applied Sciences supported this work.

Institutional Support

The authors thank the Department of Biochemistry technical staff for instrument maintenance.

Contributions Not Qualifying for Authorship

Dr. Ricardo Sousa provided guidance on spectral deconvolution methods.

References

1. Kelly SM, Jess TJ, Price NC. How to study proteins by circular dichroism. *Biochimica et Biophysica Acta*. 2005;1751(2):119-139.
2. Greenfield NJ. Using circular dichroism spectra to estimate protein secondary structure. *Nature Protocols*. 2006;1(6):2876-2890.
3. Joshi V, Shivach T, Kumar V, Yadav N, Rathore A. Circular dichroism spectroscopy as a tool for monitoring aggregation in monoclonal antibody therapeutics. *Analytical Chemistry*. 2014;86(23):11606-11613.
4. Pace CN, Shaw KL. Linear extrapolation method of analyzing solvent denaturation curves. *Proteins: Structure, Function, and Bioinformatics*. 2000;41(S4):S72-S85.
5. Lakowicz JR. *Principles of Fluorescence Spectroscopy*. 3rd ed. New York: Springer; 2006.
6. Sreerama N, Woody RW. Estimation of protein secondary structure from circular dichroism spectra. *Analytical Biochemistry*. 2000;287(2):252-260.
7. Moriyama Y, Ohta D, Hachiya K, Mitsui Y, Takeda K. Fluorescence behavior of tryptophan residues of bovine serum albumin. *Journal of Protein Chemistry*. 1996;15(3):265-272.
8. Michnik A. Thermal stability of bovine serum albumin: differential scanning calorimetry study. *Journal of Thermal Analysis and Calorimetry*. 2003;71(2):509-519.
9. Privalov PL, Dragan AI. Microcalorimetry of biological macromolecules. *Biophysical Chemistry*. 2007;126(1-3):16-24.
10. Wen J, Arthur K, Engel R, Moeller T. Introduction to biophysics for the development of biopharmaceuticals. *Journal of Pharmaceutical Sciences*. 2012;101(3):955-964.
11. Woody RW. Circular dichroism. *Methods in Enzymology*. 1995;246:34-71.
12. Sreerama N, Venyaminov SY, Woody RW. Estimation of number of secondary structure components from protein CD spectra. *Analytical Biochemistry*. 1999;273(1):50-58.
13. Johnson WC. Protein secondary structure and circular dichroism. *Proteins: Structure, Function, and Bioinformatics*. 1990;7(3):205-214.
14. Whitmore L, Wallace BA. DICHROWEB: an online circular dichroism analysis program. *Bioinformatics*. 2004;20(15):2459-2460.