


International Journal of Advanced Biochemistry Research



ISSN Print: 2617-4693
ISSN Online: 2617-4707
NAAS Rating (2026): 5.29
IJABR 2026; 10(9): 440-443
www.biochemjournal.com
Received: 17-07-2026
Accepted: 19-08-2026

Dr. Pankaj Patidar
M.V.Sc. Student, Batch 2025,
Department of Veterinary Surgery,
Co.V.Sc and A.H., Rewa,
N.D.V.S.U., Madhya Pradesh,
India

Dr. Mansi Shukla
Associate Professor, Department of
Veterinary Anatomy, Co.V.Sc and
A.H., Rewa, N.D.V.S.U., Madhya
Pradesh, India

Dr. Rajesh Ranjan
Assistant Professor, Department of
Veterinary Anatomy, Co.V.Sc and
A.H., Rewa, N.D.V.S.U., Madhya
Pradesh, India

Dr. Shivam Gupta
M.V.Sc. Student, Batch 2025,
Department of Veterinary Surgery,
Co.V.Sc and A.H., Rewa,
N.D.V.S.U., Madhya Pradesh,
India

Dr. Patil Shiva Kumar
M.V.Sc. Student, Batch 2025,
Department of Veterinary Surgery,
Co.V.Sc and A.H., Rewa,
N.D.V.S.U., Madhya Pradesh,
India

Dr. Rahul Mewada
M.V.Sc. Student, Batch 2025,
Department of Veterinary Surgery,
Co.V.Sc and A.H., Rewa,
N.D.V.S.U., Madhya Pradesh,
India

Dr. Mokkarla Amarnath Reddy
M.V.Sc. Student, Batch 2025,
Department of Veterinary Surgery,
Co.V.Sc and A.H., Rewa,
N.D.V.S.U., Madhya Pradesh,
India

Dr. Sathya Lakshmi K
M.V.Sc. Student, Batch 2025,
Department of Veterinary Surgery,
Co.V.Sc and A.H., Rewa,
N.D.V.S.U., Madhya Pradesh,
India

Corresponding Author:
Dr. Mansi Shukla
Associate Professor, Department of
Veterinary Anatomy, Co.V.Sc and
A.H., Rewa, N.D.V.S.U., Madhya
Pradesh, India

Alizarin Red S. technique and its application in veterinary sciences: A review

Pankaj Patidar, Mansi Shukla, Rajesh Ranjan, Shivam Gupta, Patil Shiva Kumar, Rahul Mewada, Mokkarla Amarnath Reddy and Sathya Lakshmi K

DOI: <https://www.doi.org/10.33545/26174693.2026.v10.i9f.10018>

Abstract

Alizarin Red S (ARS) is a primary calcium-binding dye utilized in veterinary sciences to visualize mineralized bone and calcium deposition in whole-mount skeletal preparations, histological tissues, and osteogenic cell cultures. This review outlines the fundamental chemical principles, laboratory workflow, and technical protocols associated with ARS staining, including specimen preparation, staining, washing, clearing, and mounting. By highlighting primary diaphysis and secondary epiphysis ossification sites across various developmental stages, ARS serves as an essential bridge connecting veterinary anatomy, embryology, radiology, and orthopedics. Understanding these normal skeletal mineralization patterns enables clinicians and researchers to accurately distinguish normal developmental structures from fractures, physal injuries, or congenital skeletal anomalies in immature animals. Furthermore, the review addresses key visual observations, methodological limitations, common technical artifacts, and strategies for protocol optimization in developmental research and surgical education.

Keywords: Skeletal mineralization, ossification centers, whole-mount staining, veterinary anatomy, developmental orthopaedics, veterinary radiology

Introduction

Alizarin Red S (ARS) is a calcium-binding dye used to demonstrate mineralized calcium deposits. In skeletal preparations it can make mineralized bone conspicuous and is therefore useful for studying the appearance and distribution of ossified elements during development. Classic skeletal methods have used ARS to stain ossified centres in embryos and to prepare cleared specimens for demonstration and photography [2, 5]. The technique connects developmental anatomy with radiology, orthopaedics and surgery. It is especially useful for understanding the normal pattern of skeletal mineralization in immature animals and for distinguishing developmental structures from lesions when correlated with anatomy and imaging [11]. A veterinary student should be able to explain what ARS demonstrates, describe the laboratory workflow, identify probable ossification sites and discuss their relevance to immature-animal orthopaedic and surgical assessment [2, 5, 11].

Principle and types

Principle

ARS binds calcium and produces red/orange-red coloration in calcium-rich mineralized deposits. In osteogenic cell cultures it demonstrates extracellular calcium deposition; in whole-mount skeletons it highlights mineralized bone after appropriate preparation and clearing. Common calcium-staining protocols use an acidic solution around pH 4.1–4.3, often pH 4.2, but exact conditions depend on specimen and protocol [7].

Table 1: Types / Applications

Type	Material	Purpose
Whole-mount skeletal staining	Embryo/fetus/small skeleton	Demonstrate mineralized skeletal pattern
Section/tissue staining	Histological tissue	Localize calcium/mineralized areas
In-vitro mineralization assay	Osteogenic cell culture	Assess calcium deposition
Developmental skeletal mapping	Multiple developmental stages	Compare appearance of ossified elements

Important Interpretation

ARS positivity indicates mineralization/calcium deposition; it is not a complete description of mature bone. A stained focus must be anatomically identified before being named a particular ossification centre.

ARS may be combined with Alcian Blue in whole-mount studies to show cartilage and mineralized bone as complementary components.

Staining intensity alone is not a direct measure of bone strength.




Fig 1: General workflow for whole-mount skeletal demonstration.

Ossification sites: Anatomical basis

An ossification center is a region in which bone formation/mineralization begins during skeletal development. Primary centres are generally associated with

the developing shaft of a long bone, whereas secondary centres occur in epiphyses or other secondary locations. Timing varies by species and bone.




Fig 2: Conceptual long-bone pattern of primary and secondary ossification centres.

Demonstration of ossification sites

1. Prepare the specimen for stain penetration.
2. ARS binds calcium in mineralized tissue.
3. Mineralized regions become red/orange-red.
4. Map each stained focus using anatomical landmarks.
5. Compare developmental stages or species when appropriate.
6. Confirm uncertain sites with anatomy, radiography or histology.

Progressive appearance and expansion of stained skeletal elements generally reflects increasing mineralization during

development of the animal. The pattern is best interpreted as a developmental map, not as a simple measurement of bone maturity.

Applications in veterinary science

Relation to veterinary anatomy

- Visualizes three-dimensional distribution of mineralized skeletal elements.
- Helps students understand embryology and skeletal development.
- Provides a teaching bridge between anatomy and radiology.

Relation to Veterinary Radiology

- Knowledge of normal ossification helps distinguish developmental mineralization from pathology.
- Supports interpretation of fetal and neonatal skeletal radiographs.
- Stained specimens can be compared with radiographic images.

Relation to veterinary surgery

The stain is primarily a laboratory/teaching/research method, not a routine intraoperative treatment. Its surgical importance comes from understanding the developmental anatomy of immature bones, particularly growth regions and bone ends.

Table 2: Importance of the technique

Surgical/orthopedic context	Why ossification knowledge matters
Fracture near bone end	Distinguishes developmental structures from fracture lines with imaging/anatomy.
Physeal/growth-region injury	Recognizes the region because growth disturbance may occur.
Congenital skeletal abnormality	Helps to understand whether an abnormal pattern may be developmental.
Bone surgery in young animals	Helps to plan approaches with awareness of developing structures.
Experimental orthopaedics	ARS may provide a research endpoint for mineralization.
Teaching surgical anatomy	Makes growth-related anatomy easier to visualize.

Developmental Research

Whole-mount staining can be used in developmental and toxicological studies to assess skeletal mineralization

patterns. Delayed, reduced or abnormal mineralization should be interpreted using controls and complementary methods rather than stain appearance alone ^[5, 11].

Observation, advantages, limitations and artefacts




Figure 3. Conceptual progression of mineralized skeletal appearance across development.

Table 3: Observations and Possible Interpretations

Observations	Possible interpretation	Examine
Red/orange-red focus	Mineralized element	Anatomical identity and age/stage
Stained shaft	Mineralized diaphysis	Expected primary centre
Stained epiphysis	Secondary mineralization	Species/age anatomy
Weak/patchy stain	Limited mineralization or technical issue	pH, penetration, fixation, stain quality
High background	Excess unbound dye	Washing/protocol
Asymmetry	True variation or artefact	Compare anatomy/replicates

Advantages

- Clear visual demonstration of mineralized skeletal elements.
- Useful for whole-mount developmental anatomy.
- Photographic documentation is possible.
- Useful for developmental comparisons.
- Links anatomy, embryology, radiology and orthopaedics.

Limitations

- Shows mineralization rather than every component of bone.
- Intensity is affected by preparation and protocol.
- Whole-mount preparation can be technically demanding.
- Poor clearing can obscure small centres.
- Species-specific developmental timing limits direct cross-species assumptions.

Table 4: Common Technical Problems

Problem	Possible cause	Approach
Poor staining	pH/penetration/protocol issue	Verify validated method and pH
High background	Excess dye/inadequate washing	Improve washing
Uneven staining	Poor penetration	Improve preparation
Loss of morphology	Over-clearing/rough handling	Reduce damage
Difficult interpretation	No developmental reference	Use species-specific anatomy

Integration with veterinary surgery

A growing animal is not simply a smaller adult. Developing bones contain growth regions and changing patterns of mineralization. Recognition of normal development is therefore important when assessing trauma, congenital defects and orthopaedic disease in young animals ^[11-15].

Table 5: Interpretation and Integration

Area	Clinical question	Role of developmental knowledge
Fracture near bone end	Fracture or normal developmental structure?	Use age-specific anatomy and imaging.
Growth plate injury	Could growth be disturbed?	Recognize the physis/growth region.
Congenital deformity	Is it developmental?	Know normal skeletal sequence.
Surgery in young animal	Could a developing structure be affected?	Plan with developmental anatomy.
Bone-healing research	Is mineralization changing?	ARS may be a validated experimental endpoint.

Conclusion

Alizarin Red S is an important calcium-binding stain for demonstrating mineralized tissue. In whole-mount skeletal preparations it provides a useful visual method for demonstrating mineralized/ossified skeletal elements during development. For veterinary students, its value extends into radiology, orthopaedics and surgery because understanding normal ossification patterns is essential when evaluating immature animals. Correct interpretation requires species, age and anatomical correlation.

References

1. Hood RCW, Neill WM. A modification of Alizarin Red S technic for demonstrating bone formation: methods of mounting, photographing and demonstrating the specimens. *Stain Technology*. 1948;23(4):209–218. DOI:10.3109/10520294809106251.

2. Gregory CA, Gunn WG, Peister A, Prockop DJ. An Alizarin Red-based assay of calcium accumulation by osteogenic cells. *Methods in Molecular Biology*. 2004;240:217–224.

3. Hall LW, Clarke KW, Trim CM. *Veterinary Anaesthesia*. 10th ed. W.B. Saunders; 2001.

4. Fossum TW, editor. *Small Animal Surgery*. 5th ed. Elsevier; 2019.

5. Dyce KM, Sack WO, Wensing CJG. *Textbook of Veterinary Anatomy*. 5th ed. Elsevier; 2017.

6. Fubini SL, Ducharme NG, editors. *Farm Animal Surgery*. 2nd ed. Elsevier; 2016.

7. Bernar A, Gebetsberger JV, Bauer M, Streif W. Optimization of the Alizarin Red S assay by enhancing mineralization of osteoblasts. *International Journal of Molecular Sciences*. 2023;24(1):723. DOI:10.3390/ijms24010723.

8. Qiu D, *et al*. *In vitro* determination of osteo-adipogenic lineage choice of bone marrow stromal/stem cells (BMSCs). *MethodsX*. 2024;12:102637. DOI:10.1016/j.mex.2024.102637.

9. Onoki T, *et al*. Modulation of osteoblastogenesis by NRF2: NRF2 activation suppresses osteogenic differentiation and enhances mineralization in human bone marrow-derived mesenchymal stromal cells. *FASEB Journal*. 2024. DOI:10.1096/fj.202400602R.

10. Huang CY, Lin WL. High-throughput quantification of osteogenesis in WJ-MSCs using Alizarin Red S, von Kossa, and Picro-Sirius red. *Future Science OA*. 2025;11(1):2567169. DOI:10.1080/20565623.2025.2567169.

11. Kim JK, Park D, Chang MC. Assessment of bone age based on hand radiographs using regression-based multi-modal deep learning. *Life*. 2024;14(6):774. DOI:10.3390/life14060774.

12. Van den Broeck M, Chen Y, Cornillie P. Age estimation of puppies based on the radiographically assessed development of ossification centres in the carpal and metacarpal regions. *Veterinary Record*. 2023;192(9). DOI:10.1002/vetr.2582.

13. Rao SD, Rajendranath N, Chandrasekhara Rao TS, Ramayya PJ, Suresh Kumar RV. Prenatal development of bones of hip, thigh and leg regions of buffalo. *Indian Journal of Veterinary Anatomy*. 2012;24(1).

14. Nissar S, Basha SH, Ramesh G, Mayakkannan T, Vimal R. Prenatal study on ossification and growth of long bones in sheep foetus. *Indian Journal of Veterinary Anatomy*. 2018;29(1).

15. Karunasri V. Sequential changes in progression of ossification in middle ear ossicles of Nellore sheep fetus. *Indian Journal of Veterinary Anatomy*. 2018;30(1).